

Premium Foley Catheters Promotion

Offer valid in the USA and Canada January 22, 2018 – June 30, 2018

Buy one Premium Foley Catheter, receive one free!

Premium Foley Catheters**

- Closed-end silicone catheter with distal balloon cuff used for intermittent drainage of urine from the bladder or contrast studies of urogenital tract
- Balloon aids in retaining the catheter inside the bladder
- Recommended to be used with wire guide stiffening stylets
- Sterile
- Radiopaque
- Atraumatic
- Autoclavable
- Color coded for easy size reference

** Use only sterile water to inflate balloon



Product Code	French Size	Length	Balloon Size	Color
FC530	5	30 cm or 12"	2 cc	Pink
FC555	5	55 cm or 22"	2 cc	Pink
FC1230	12	30 cm or 12"	10 cc	White

For more information, contact your Distribution Partner or SurgiVet® Sales Account Manager