

Quarterly Promotion

Offer valid in the USA and Canada April 1, 2018 – June 30, 2018

Advisor® Vital Signs Monitor Offer

Purchase an Advisor® Vital Signs Monitor* and receive a complimentary Roll Stand (3387) – a \$374 value

- Monitor is available in multiple configurations. It includes three monitoring modes that adjust ECG sensitivity, cuff inflation pressure and alarm limits to tailor to the patient
 - Standard parameters include ECG, digital pulse oximetry, heart rate, non-invasive blood pressure and a built-in battery
 - Optional parameters include end tidal CO₂, respiratory rate, 2-channel temperature, 2-channel invasive blood pressure, and an integrated printer
- The roll stand has two locking casters and a mounting plate with locking device for the easy attachment or removal of the monitor. The utility basket provides convenient storage. Dimensions: 45" H x 23" W.

* The following Advisor product codes are included in the offer: V9200, V9200P, V9200IBP/TEMP, V9201, V9201P, V9203, V9204



V9203



3387


Refurbished Medfusion® 3500 Syringe Pump Offer

Receive 25% off the purchase a refurbished Medfusion® 3500 syringe pump (3500-0600-82) and receive a complimentary case of 25 Medex® Mangum™ Micro Fluid Delivery sets (MX448FMG) – a \$159 value

- The syringe pump is an intuitive system designed to enhance safety while maintaining simplicity. Precise delivery of fluids and medications allows for the implementation of a wider range of treatment modalities.
- The Mangum™ micro delivery sets are designed to enable primary IV administration of metered infusion from a Medfusion® syringe pump.



3500-0600-82



MX448FMG

HOTLINE® Fluid Warmer Offer

Purchase a HOTLINE® Fluid Warmer (V707000) and receive a complimentary case of 30 blood and IV fluid warming sets with injection port (L-70) – a \$243 value

- Allows fluid or blood to be warmed to the appropriate range of 95-105° Fahrenheit
- A patented patient line with triple lumen tubing maintains a layer of 107.6° Fahrenheit solution circulating around IV line all the way to the patient
- Warm and ready to use in four minutes (from a cold start)
- Accommodates a slow drip to a bolus
- Can be mounted to an IV pole
- Remanufactured: one-year warranty
- Dimensions: 9.5" H x 8.3" W x 7.0" D
- Weight: 7.95 lbs.



V707000



L-70

**For more information, contact your Distribution Partner or
SurgiVet® Sales Account Manager**